

DNA Technical Information

Adapted by Lori Yang from www.biosyn.com/t_dna.htm

I. Calculations

A. Using Absorbance at 260 nm to Measure DNA Concentration

A useful estimation of DNA concentration relates to the amount of oligonucleotide which, when dissolved in 1 mL of water, results in an absorbance of 1 when measured at 260 nm in a 1 cm path length cuvette. This is often simply called the A_{260} of a sample. The actual concentration can range from 39 $\mu\text{g/mL}$ (for a homopolymer of C) to 20 $\mu\text{g/mL}$ (for a homopolymer of A). For most practical experiments, an A_{260} of 1.0 represents approximately 33 mg of oligo with an equal mixture of the four bases.

A_{260} conversion factors

$A_{260} = 1.0 \rightarrow 33 \mu\text{g/mL ssDNA}$
 $\rightarrow 40 \mu\text{g/mL ssRNA}$
 $\rightarrow 50 \mu\text{g/mL dsDNA}$

Molecular Weight of an Oligonucleotide

$MW_{\text{oligo}} = ((A \times 312.2) + (G \times 328.2) + (C \times 288.2) + (T \times 303.2)) - 61$

where A,C,G,T represent the number of A's, C's, G's and T's in an oligo.

Molar conversions

1 μg of 1,000 bp DNA = 1.52 pmol (3.03 pmoles of ends)

1 μg of pBR322 DNA = 0.36 pmol DNA

1 pmol of 1,000 bp DNA = 0.66 μg

1 pmol of pBR322 DNA = 2.78 μg

Melting Temperature (Nucleic Acid Hybridization)

Up to 25 bp:

$$T_m = 4^\circ\text{C} (G+C) + 2^\circ\text{C} (A+T)$$

More than 25 bp:

$$T_m = 81.5^\circ\text{C} + 16.6 \log M + 0.41 \% (G+C) - 500 / n - 0.61 (\% \text{ formamide})$$

$M = [\text{Na}^+]$ in moles/liter ; n = length of shortest chain in duplex

B. Other Information

Resuspension Buffers

1. Sterile Water (dd H_2O)
2. TE Buffer (10mM Tris-HCl, 1mM EDTA) pH 7.5

DNA Storage Conditions and Stability

Lyophilized (-20°C) = 6 months to several years

Lyophilized (25°C) = 2 months to 1 year

Dissolved (-20°C) = 1 month to 6 months

Dissolved (25°C) = 1 week to 3 months

DNA Conformations

Helix type	Direction of rotation	Residues per turn	Rotation per residue	Helix rise per residue	Helix pitch
A	Right	11	33°	2.55 Å	28 Å
B	Right	10	36°	3.4 Å	34 Å
Z	Left	12	-30°	3.7 Å	45 Å

PAGE Purification Information

% Gel	Dye migration on a denaturing gel (7 M UREA), 1x TBE		
	Separation range (bases)	Bromophenol Blue	Xylene Cyanol
5.0	100-200	35	140
8.0	80-150	20	75
12.0	40-100	12	53
15.0	12-80	10	44
20.0	8-60	8	24

References

1. Suggs, S., et. al. Proc. Natl. Acad. Sci 78:6613('81)
2. George H. Keller, Mark M. Manak, DNA probes; p 15; M Stockton Press, '89.